


O.R. COUPLERS

Repair Parts Lists

Allied
HEALTHCARE PRODUCTS INC.

CHEMETRON Medical Division

ENGR. # S171-122-001
FORM # 11-00-2001
DATE JAN. 2020 REV. E



For Valve Body O-Ring Service Use Chemetron O-Ring Service Tool 64-90-2060

Refer to Form 64-00-2008 Medical Gas Outlets for Maintenance Procedures

For Disassembly of Coupler Valve Body and Poppet Replacement Use Repair Tool 65-10-1001 Form Number 64-00-1010

NOTE: Lubricate all O-Rings with 64-90-2111 before installation.

Allied
HEALTHCARE PRODUCTS INC.
GLOBAL SUPPORT OF LIFE™

CHEMETRON Medical Division

1720 Sublette Avenue
St. Louis, MO 63110
Telephone: (314) 771-2400
Toll Free: (800) 444-3954

INSTRUMENTS OF CARE